

Laser Hair Removal FAQs

Q? What does Matrix use for Laser Hair Removal?

A: We use Palomar Stralux 500 that is medical grade laser system. It is an extremely well build versatile platform with large range of applications. This USA made system is the Rolls-Royce of lasers. This is the best laser system in the world build to perfection.

Q? How many treatments are necessary?

A: Some clients will benefit from only 2 to 3 treatment sessions. Sometimes more treatments are necessary to achieve permanency. Pulsed light energy has its greatest effect on growing, or anagen hair follicles; as the dormant, or telogen hairs become active, usually within 3 weeks to 3 months, some new hair growth will be observed.

Q? What type of hair & skin can be treated?

A: Specific wavelengths of light can effectively target both finer and coarser hair as well as lighter hair, on all skin types.